Product Instructions

- ⓔⓃ Neogen® Molecular Detection – Right Now™ Matrix Control
- ⓕⓇ Contrôle de matrice Neogen® Molecular Detection – Right Now™
- ⓓⓔ Neogen® Right Now™ Matrixkontrolle für molekulare Detektion
- ⓔⓈ Control de matriz Right Now™ para detección molecular de Neogen®
- ⓅⓉ Detecção Molecular Neogen® – Controle de matriz Right Now™
- ⓙⓐ Neogen® 分子検出-現在™マトリックス制御
- Ⓩⓗ Neogen® 分子检测 – Right Now™ 基质对照
- ⓉⓌ Neogen® 分子檢測 – Right Now™ 基質對照
- Ⓣⓓ ตัวควบคุมคุณภาพของเมทริกซ์ (Matrix Control) ของการตรวจหาระดับโมเลกุลของ Neogen® – Right Now™
- Ⓚⓐ Neogen® 분자 검출 – Right Now™ Matrix Control

Product Instructions

Neogen Molecular Detection – Right Now Matrix Control

Product Description and Intended Use

The Neogen Molecular Detection – Right Now Matrix Control (RNMC96) is intended for use in pair with the Neogen® Molecular Detection Assay, *Listeria* Right Now™ (MDALRN96) and the Neogen® Molecular Detection System™ (MDS100). Its purpose is to help users demonstrate the robustness of the MDALRN96 assay in their specific processing environment during method implementation.

The Neogen Molecular Detection – Right Now Matrix Control employs loop-mediated isothermal amplification (LAMP) to rapidly amplify nucleic acid sequences with high specificity and sensitivity, combined with bioluminescence to detect the amplification. Valid results are reported in real time while inhibited results will be displayed after the assay is completed.

It is generally accepted that the sample matrix may interfere with any test method. To help customers evaluate the suitability of the Neogen Molecular Detection Assay – *Listeria* Right Now under different environmental conditions, Neogen has developed the Molecular Detection – Right Now Matrix Control kit. Use the Right Now Matrix Control to verify that the sample matrix does not interfere with the assay. It is recommended to test representative samples (i.e. samples from different surface types, equipment or if applicable, different sanitation procedures) during method adoption or whenever cleaning and sanitation protocols have changed.

The Neogen Molecular Detection Assay – *Listeria* Right Now (MDALRN96) and Neogen Molecular Detection – Right Now Matrix Control (RNMC96) have only been evaluated for use with environmental samples from visibly clean surfaces, free of excess sanitizers (post-cleaning/sanitation and pre-production).

The Neogen Molecular Detection – Right Now Matrix Control has not been evaluated with all possible food processing environments. Neogen has not documented the use of this product in industries other than food or beverage. Neogen recommends evaluation of the method in the user's environment using a sufficient number of samples for particular environmental challenges to ensure that the method meets the user's criteria.

The Neogen Molecular Detection – Right Now Matrix Control is intended for use in a laboratory environment by professionals trained in laboratory techniques.

The Neogen® Molecular Detection Instrument is intended for use with samples that have undergone heat treatment during the assay lysis step, which is designed to lyse the cells and release the nucleic acids from the organisms present in the sample. Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should not be inserted into the Neogen Molecular Detection Instrument.

Neogen Food Safety is certified to International Standards Organization (ISO) 9001 for design and manufacturing.

Neogen Molecular Detection – Right Now Matrix Control contains 96 tests described in Table 1.

Table 1. Kit Components

Item	Identification	Quantity	Contents	Comments
Right Now Matrix Control (RNMC) Reagent Tubes	Clear tubes	96 (12 strips of 8 tubes)	Lyophilized nucleic acids fragment, specific amplification and detection mix	Ready-to-use
Extra caps	Clear caps	96 (12 strips of 8 caps)		Ready-to-use

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Molecular Detection System and the Neogen Molecular Detection – Right Now Matrix Control. Retain the safety instructions for future reference.

⚠ **WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ **CAUTION:** Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

⚠ **WARNING**

Do not use the Neogen Molecular Detection – Right Now Matrix Control in the diagnosis of conditions in humans or animals.

The user must train personnel in current proper testing techniques, for example: Good Laboratory Practices, ISO/IEC 17025(1), or ISO 7218(2).

To reduce the risks associated with a false-negative result leading to the release of contaminated surfaces:

- Follow the protocol and perform the tests exactly as stated in the product instructions.
- Always use Good Laboratory Practices to transfer samples and reagents during execution of the method.
- Store the Neogen Molecular Detection – Right Now Matrix Control as indicated on the package and in the product instructions.
- Always use the Neogen Molecular Detection – Right Now Matrix Control by the expiration date.
- Use the Neogen Molecular Detection – Right Now Matrix Control to assess sample interference during method adoption or when cleaning and sanitation protocols have changed.
- Use the Neogen Molecular Detection – Right Now Matrix Control only with surfaces, sanitizers, and protocols that have been validated internally or by a third party. The Neogen Molecular Detection – Right Now Matrix Control is not designed for testing soiled surfaces or for food matrix testing.

To reduce the risks associated with exposure to chemicals and biohazards:

- Perform pathogen testing in a properly equipped laboratory under the control of trained personnel.
- Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection, while handling reagents and contaminated samples.
- Avoid contact with the contents of the sample, lysis tubes, and reagent tubes after amplification.
- Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should not be inserted into the Neogen Molecular Detection Instrument.

To reduce the risks associated with cross-contamination while preparing the assay:

- Always wear gloves (to protect the user and prevent introduction of nucleases).

To reduce the risks associated with environmental contamination:

- Follow current industry standards for disposal of contaminated waste.

⚠ **CAUTION**

To reduce the risks associated with exposure to hot liquids:

- Do not exceed the recommended temperature setting on the heater.
- Do not exceed the recommended heating time.
- Use an appropriate, calibrated thermometer to verify the Neogen® Molecular Detection Heat Block Insert temperature (e.g., a partial immersion thermometer or digital thermocouple thermometer, not a total immersion thermometer). The thermometer must be placed in the designated location in the Neogen Molecular Detection Heat Block Insert.

To reduce the risks associated with cross-contamination while preparing the assay:

- Change gloves prior to reagent pellet hydration.
- Use of long-tip, sterile, aerosol barrier (filtered), molecular biology grade pipette tips are recommended.
- Use a new pipette tip for each sample transfer.

- Use Good Laboratory Practices to transfer the sample to the lysis tube. To avoid cross-contamination between samples, the user may choose to clean the shaft of the pipette with a wipe with 1–5% bleach and then with a clean paper towel or use extended length pipette tips.
- Use a molecular biology workstation containing a germicidal lamp where available. Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1–5% (v:v in water) household bleach solution or nucleic acid removal solution.

To reduce the risks associated with a false-positive result:

- Never open tubes post-amplification.
- Always dispose of the contaminated tubes by soaking in a 1–5% (v:v in water) household bleach solution for 1 hour, and away from the assay preparation area.
- Never autoclave reagent tubes post-amplification.

Consult the Safety Data Sheet for additional information and local, state, and national regulations for disposal.

If you have questions about specific applications or procedures, please visit our website at neogen.com or contact your local Neogen representative or distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at neogen.com or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, laboratory technique, and the sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet their customers' and suppliers' requirements.

As with any test method, results obtained from the use of any Neogen Food Safety product do not constitute a guarantee of the quality of the samples or processes tested.

To help customers evaluate the *Listeria* Right Now assay for various environmental conditions, Neogen has developed the Neogen Molecular Detection – Right Now Matrix Control kit. Use the Neogen Molecular Detection – Right Now Matrix Control (RNMC96) to determine that the matrix does not interfere with the Neogen Molecular Detection Assay – *Listeria* Right Now. Test several samples, representative of the environment matrix (i.e., samples obtained from different surface types, equipment, or if applicable different sanitation procedures), during any validation period when adopting the Neogen method or when cleaning and sanitation protocols have changed.

Limitation of Warranties / Limited Remedy

Except as expressly stated in a limited warranty section of individual product packaging, Neogen disclaims all express and implied warranties, including but not limited to, any warranties of merchantability or fitness for a particular use. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen Food Safety representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

Neogen will not be liable for any loss or damages, whether direct, indirect, special, incidental, or consequential damages, including but not limited to lost profits. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage and Disposal

Store the Neogen Molecular Detection – Right Now Matrix Control at 2–8°C. Do not freeze. Keep the kit away from light during storage.

After opening the kit, check that the foil pouch is undamaged. If the pouch is damaged, do not use.

After opening, unused reagent tubes should always be stored in the resealable pouch with the desiccant inside to maintain stability of the lyophilized reagents. Store resealed pouches at 2–8°C for no longer than 60 days.

Do not use Neogen Molecular Detection – Right Now Matrix Control past the expiration date. The expiration date and lot number are noted on the outside label of the box.

When testing is complete, follow national, state, and local regulations for the disposal of contaminated waste. Consult the Safety Data Sheet for additional information.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1–5% (v:v in water) household bleach solution or nucleic acid removal solution.

The user should complete the Neogen Molecular Detection System operator qualification (OQ) training, as described in the “Installation Qualification (IQ)/Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System” document (3).

Preparation of Equipment

1. Preparation of the Neogen® Molecular Detection Speed Loader Tray
 - 1.1 Wet a cloth or disposable towel with a 1–5% (v:v in water) household bleach solution and wipe the Neogen Molecular Detection Speed Loader Tray.
 - 1.2 Rinse the Neogen Molecular Detection Speed Loader Tray with water.
 - 1.3 Use a disposable towel to wipe the Neogen Molecular Detection Speed Loader Tray dry.
 - 1.4 Ensure the Neogen Molecular Detection Speed Loader Tray is dry before use.
 2. Preparation of the Neogen® Molecular Detection Chill Block Insert
 - 2.1 Place the Neogen Molecular Detection Chill Block Insert directly on the laboratory bench. Use the block at ambient laboratory temperature.
 3. Preparation of the Neogen® Molecular Detection Heat Block Insert
 - 3.1 Place the Neogen Molecular Detection Heat Block Insert in a dry double block heater unit. Turn on the dry block heater unit and set the temperature to allow the Neogen Molecular Detection Heat Block Insert to reach and maintain a temperature of $100 \pm 1^{\circ}\text{C}$.
- Note:** Depending on the heater unit, allow approximately 30 minutes for the Neogen Molecular Detection Heat Block Insert to reach temperature. Using an appropriate, calibrated thermometer (e.g., a partial immersion thermometer, digital thermocouple thermometer, not a total immersion thermometer) placed in the designated location, verify that the Neogen Molecular Detection Heat Block Insert is at $100 \pm 1^{\circ}\text{C}$.
4. Preparation of the Neogen Molecular Detection Instrument
 - 4.1 Launch the Neogen® Molecular Detection Software and log in. Contact your Neogen Food Safety representative to ensure you have the most updated version of the software.
 - 4.2 Turn on the Neogen Molecular Detection Instrument.
 - 4.3 Create a new run or edit a run by selecting a well for each sample. Refer to the Neogen Molecular Detection System User Manual for details.

- 4.4 For the Neogen Molecular Detection – Right Now Matrix Control (RNMC96), select the MC-RN icon on the Neogen Molecular Detection Software and assign to each well.
- 4.5 For the *Listeria* Right Now assay (MDALRN96), select the LRN icon on the Neogen Molecular Detection Software and assign to each well.
- 4.6 Use the same sample ID for both the LRN and MC-RN wells. The software will link the wells and assess the results of both reactions in parallel. The outcome of the LRN reaction will depend on the MC-RN result to determine if the sample is inhibited.

Note: The Neogen Molecular Detection Instrument must reach and maintain a temperature of 60°C before inserting the Neogen Molecular Detection Speed Loader Tray with reaction tubes. This heating step takes approximately 20 minutes and is indicated by an orange light on the instrument's status bar. When the instrument is ready to start a run, the status bar will turn green.

Sample Preparation and Lysis

For detailed instructions, please refer to the Product Instructions for the specific Right Now assay (e.g MDALRN96) being used.

Note: Extended-length pipette tips (that are longer than standard pipette tips) are recommended. In case the pipette shaft touches the Right Now Sampler tube, eject the pipette tip, then use a towel sprayed with 1–5% bleach to thoroughly wipe the shaft of the pipette to prevent cross-contamination between samples.

Amplification


Transfer sample lysate to Neogen Molecular Detection – Right Now Matrix Control tubes as described below:

1. Use one Neogen Molecular Detection – Right Now Matrix Control (RNMC96) tube for each sample.
 - 1.1 Neogen Molecular Detection – Right Now Matrix Control tube strips can be cut to desired tube number. Select the number of individual tubes or 8-tube strips needed.
 - 1.2 Place Neogen Molecular Detection –Right Now Matrix Control tubes in an empty rack.
 - 1.3 Avoid disturbing the reagent pellets from the bottom of the tubes.
2. To avoid cross-contamination, decap one Neogen Molecular Detection – Right Now Matrix Control tube strip at a time, and use a new pipette tip for each transfer step.
3. Transfer lysate to Neogen Molecular Detection – Right Now Matrix Control tubes (RNMC) as described below:
 - 3.1 Use the Neogen® Molecular Detection Cap/Decap Tool-Reagent to decap the RNMC tubes — one RNMC strip at a time. Discard cap.
 - 3.2 Transfer 50 µL of sample lysate from the upper half of the liquid in the Neogen Lysis Solution tube (avoiding any precipitate) into corresponding Neogen Molecular Detection – Right Now Matrix Control tube. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down five times to dissolve lyophilized reagent.


Note: The same sample lysate should be inoculated twice. Transfer 50 µL to both a LRN Reagent tube and an RNMC tube so the same lysate is analyzed in parallel in both tubes.

 - 3.3 Repeat step 3.2 until each individual sample lysate has been added to a corresponding Neogen Molecular Detection – Right Now Matrix Control tube in the strip.
 - 3.4 Cover the Neogen Molecular Detection – Right Now Matrix Control tubes with the provided extra cap and use the rounded side of the Neogen Molecular Detection Cap/Decap Tool-Reagent to apply pressure in a back and forth motion ensuring that the cap is tightly applied.
 - 3.5 Repeat steps 3.1 to 3.4 as needed, for the number of selected environmental sample matrices to be tested.

- Load capped tubes into a clean and decontaminated Neogen Molecular Detection Speed Loader Tray, then close and latch the lid.


Transfer 50 μ L of sample lysate from the upper $\frac{1}{2}$ of the cluster tube into the RNMC Tube. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.

Load capped tubes into a clean and decontaminated Neogen Molecular Detection Speed Loader Tray. Close and place it into the Neogen MDS. Close the lid to start the assay

- Review and confirm the configured run in the Neogen Molecular Detection Software.

- Click the start button in the software and select instrument for use. The selected instrument's lid automatically opens.
- Place the Neogen Molecular Detection Speed Loader Tray into the Neogen Molecular Detection Instrument and close the lid to start the assay. Results are provided within 45 minutes, although a valid RNMC result may be detected sooner.
- After the assay is complete, remove the Neogen Molecular Detection Speed Loader Tray from the Neogen Molecular Detection Instrument and dispose of the tubes by soaking in a 1–5% (v:v in water) household bleach solution for one hour and away from the assay preparation area.

Note: To minimize the risk of false positives due to cross-contamination, never open reagent tubes containing amplified DNA. Always dispose of sealed reagent tubes by soaking in a 1–5% (v:v in water) household bleach solution for one hour, and away from the assay preparation area.

Results and Interpretation

An algorithm interprets the light output curve resulting from the detection of the nucleic acid amplification. Results are analyzed automatically by the software and are color-coded based on the result. A Valid or Inhibited result is determined by analysis of a number of unique curve parameters. Valid Neogen Molecular Detection – Right Now Matrix Control results indicate no inhibitory effect from the matrix sample. Valid results are reported in real time while inhibited results will be displayed after the run is completed.

If during software configuration, both the LRN well and the MC-RN well are assigned the same ID and inhibition is detected, both wells will be marked as Inhibited. Note that a positive sample well overrules the inhibited call, and will be marked as positive.

If the LRN well and the MC-RN well are assigned different IDs and inhibition is detected, the LRN well may be marked as negative, while the MC-RN well will be marked as inhibited.

Note: A no-amplification sample will not give a zero reading as the system and Neogen Molecular Detection – Right Now Matrix Control amplification reagents have a “background” relative light unit (RLU) reading.

Neogen recommends the user to repeat the assay for any inhibited samples as follows:

1. For troubleshooting purposes only and to confirm the inhibition: Re-test the original sample lysate in parallel with a diluted sample.
2. To prepare the diluted sample: Transfer 20 µL of the original sample lysate from the cluster tube into a microcentrifuge tube and add 80 µL of lysis buffer reagent solution (1:5 dilution).
3. Test the samples: Transfer 50 µL of both the original and diluted sample lysate into the Neogen Molecular Detection – Right Now Matrix Control tubes to assess whether the interference is reproducible and persists.

If you have questions about specific applications or procedures, please visit our website at neogen.com or contact your local Neogen representative or distributor.

Table 2: Result symbols for samples and interpretation as provided by the Neogen Molecular Detection System Software.

Well Type	Well Result Symbol	Result	Interpretation
Right Now Matrix Control MC-RN		Valid	The RN Matrix Control is valid.
Right Now Matrix Control MC-RN		Inhibited*	The sample matrix was inhibitory for the assay. A re-test or optimization of the sample preparation may be necessary.
Right Now Matrix Control MC-RN		Error*	No bioluminescence was detected. A re-test may be required.

*Refer to the run report for additional details on the results. If needed, contact your local Neogen representative for troubleshooting assistance.

References

¹ NF EN ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

² NF EN ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.

³ Neogen Installation Qualification (IQ) Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System.

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Referencias

- ¹ NF EN ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
- ² NF EN ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
- ³ Neogen Installation Qualification (IQ) Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System.

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Referências

- ¹ NF EN ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
- ² NF EN ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
- ³ Neogen Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System.

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