



**3M** Science.  
Applied to Life.™

# Free up your focus.

## **3M™ Molecular Detection Assay 2 - *Campylobacter* and 3M™ *Campylobacter* Enrichment Broth**

**Primary detection method used by the USDA Food Safety and Inspection Service<sup>1</sup>**

There's now a better way to test for bacteria associated with poultry production. Find new efficiencies and simplify *Campylobacter* testing with 3M™ Molecular Detection Assay 2 - *Campylobacter*. Featuring all the benefits of the 3M™ Molecular Detection System, this assay and one-step enrichment broth are easy to use and accurate. It's significantly faster than alternatives like PCR, immunoassay and culture methods.

- Efficiently monitors for *Campylobacter*
- Simplified enrichment prep and protocol
- Easy to use and train
- Streamlined workflow
- Reliable
- Next-day results



## Proven performance.

**3M™ Molecular Detection Assay 2 - *Campylobacter***  
 While *Campylobacter* may not be your most frequent poultry test, it can't be ignored either — this pathogen causes an estimated 1.5 million illnesses each year in the U.S.<sup>2</sup> Now you can test for *Campylobacter* with a solution from 3M.

The 3M™ Molecular Detection System with 3M™ Molecular Detection Assay 2 - *Campylobacter* and 3M™ Molecular Detection Assay 2 - *Salmonella* is a complete solution that allows you to monitor poultry for both bacteria simultaneously. A single post-enrichment protocol for these pathogens means you can perform up to 96 tests of multiple types in one 60-minute run. And when compared using metrics such as time to result, throughput, unit size and workflow, 3M Molecular Detection Assay 2 - *Campylobacter* outperforms other methods and manufacturers.

### 3M™ *Campylobacter* Enrichment Broth

- Easy to use
- Simplified enrichment — just add sterile water
- No supplements, microaerophilic incubation, organic solvents or autoclaving the broth




Get results you can trust at  
[3M.com/foodsafety/MDA2](http://3M.com/foodsafety/MDA2).



**3M Food Safety**  
 3M Center, Building 275-5W-05  
 St. Paul, MN 55144-1000 U.S.A.

Phone 1-800 328-6553  
 Web [3M.com/foodsafety](http://3M.com/foodsafety)

## Simple enrichment and streamlined workflow.



## Ordering Information

Description	Cat. No.	Contents
3M™ Molecular Detection Assay 2 - <i>Campylobacter</i>	MDA2CAM96	96 tests/kit
3M™ <i>Campylobacter</i> Enrichment Broth	CE250	1 bottle
3M™ Molecular Detection Assay 2 - <i>Salmonella</i>	MDA2SAL96	96 tests/kit
3M™ Bird Rinse with 400 mL Buffered Peptone Water Broth	BRBPW400	35/box
3M™ Bird Rinse Plain Sample Bag	BRB3500	500/box

<sup>1</sup>United States Department of Agriculture Food Safety and Inspection Service (2021). FSIS Constituent Update, 24 (31). Retrieved from [ConstUpdate\\_04.02.2021\\_.pdf \(usda.gov\)](https://www.fsis.usda.gov/constituent-update/04-02-2021).

<sup>2</sup>US Centers for Disease Control and Protocol Yellow Book, Ch. 4 (2020)